

Textbook of Echocardiography for Intensivists and Emergency Physicians

Armando Sarti
F. Luca Lorini
Editors

Second Edition

 Springer

Textbook of Echocardiography for Intensivists and Emergency Physicians

Armando Sarti • F. Luca Lorini
Editors

Textbook of Echocardiography for Intensivists and Emergency Physicians

Second Edition

 Springer

Editors

Armando Sarti
Rianimazione
Ospedale Santa Maria Nuova
Firenze
Italy

F. Luca Lorini
Dipartimento di Anestesia e
Rianimazione
Az. Osp. Ospedali Riuniti di Bergamo
Bergamo
Italy

ISBN 978-3-319-99890-9 ISBN 978-3-319-99891-6 (eBook)
<https://doi.org/10.1007/978-3-319-99891-6>

© Springer Nature Switzerland AG 2012, 2019

This work is subject to copyright. All rights are reserved by the Publisher, whether the whole or part of the material is concerned, specifically the rights of translation, reprinting, reuse of illustrations, recitation, broadcasting, reproduction on microfilms or in any other physical way, and transmission or information storage and retrieval, electronic adaptation, computer software, or by similar or dissimilar methodology now known or hereafter developed.

The use of general descriptive names, registered names, trademarks, service marks, etc. in this publication does not imply, even in the absence of a specific statement, that such names are exempt from the relevant protective laws and regulations and therefore free for general use.

The publisher, the authors, and the editors are safe to assume that the advice and information in this book are believed to be true and accurate at the date of publication. Neither the publisher nor the authors or the editors give a warranty, expressed or implied, with respect to the material contained herein or for any errors or omissions that may have been made. The publisher remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.

This Springer imprint is published by the registered company Springer Nature Switzerland AG
The registered company address is: Gewerbestrasse 11, 6330 Cham, Switzerland

Preface

The first edition of this book, *Echocardiography for Intensivists*, was published six years ago. At that time, the use of ultrasound machines was not yet a routine for many intensivists around the world. Today, many doctors, who are involved in the management of unstable or critically ill patients, are simply unable to practice without it. Echocardiography is now firmly considered to be the most useful diagnostic tool for the evaluation of patients with acute cardiovascular disorders in the ICU, the operating theatre and the emergency department. Ultrasound examination of the heart and great vessels provides reliable functional anatomy and hemodynamic assessment of the cardiovascular system. There is no alternative way to assist and guide the clinicians so that, within just a few minutes, they may understand the cause of cardiocirculatory shock and accordingly treat hemodynamic derangement at the bedside.

Today, more than 60 years after its initial conception, echocardiography is established enough to be in a widespread use even outside the cardiology unit. Ultrasound assessment of the cardiovascular system is continuously being refined according to advances in research and technology. However, the appropriate clinical application of this unique diagnostic technique in the acute setting requires specific education with constant updates and training.

The first edition of the book was developed from awareness of the need for a text specifically written for intensivists, emergency physicians and anaesthesiologist who wished to incorporate the ultrasound technique into their clinical practice. Thanks to the contributions made by all the authors, the first edition of the book has been well-received all around the world, and a Chinese language edition has also been published. This second edition has been completely revised and expanded with the aim of becoming a *Textbook of Echocardiography for Intensivist and Emergency Physicians*. Each existing chapter has been updated, with the addition of many new figures and references, and ten entirely new chapters have been added.

We wish to express our gratitude to all of the authors of both editions for their invaluable contributions, and we would also like to thank Springer for their skilful production of this textbook. Last, but not least, we are indebted to our families for their understanding and support.

Finally, we wish to dedicate our efforts to all the patients who have to face critical medical conditions and to the doctors, nurses and caregivers who do their best to provide them with the care they need.

We will feel greatly rewarded if this book helps to improve, in some way, the care of such patients.

Firenze, Italy
Bergamo, Italy

Armando Sarti
F. Luca Lorini

Contents

Part I Ultrasound and Use of the Echo Machine

- 1 Essential Physics of Ultrasound and Use of the Ultrasound Machine** 3
Dionisio F. Colella, Paolo Prati, and Armando Sarti

Part II Standard Echocardiographic Examination

- 2 Ultrasound Morphology of the Heart: Transthoracic Examination** 21
Armando Sarti, Simone Cipani, and Costanza Innocenti
- 3 Transthoracic Echocardiography in the ICU: The Patient Who Is Difficult to Study** 37
Piercarlo Ballo
- 4 Ultrasound Morphology of the Heart: Transesophageal Examination.** 47
Ferdinando Luca Lorini, Carlo Sorbara, and Sergio Cattaneo
- 5 Three-Dimensional Echocardiography** 57
Mauro Pepi and Gloria Tamborini
- 6 Speckle Tracking Echocardiography (STE) for Left and Right Ventricles** 71
Simone Cipani, Anna Viappiani, and Armando Sarti

Part III Essential Functional Echo-Anatomy

- 7 The Left Ventricle: Frequent Imaging in ICU Patients.** 81
Simone Cipani, Claudio Poli, and Silvia Marchiani
- 8 The Right Ventricle and Pulmonary Artery: Frequent Imaging in ICU Patients** 101
Luigi Tritapepe, Mario Mezzapesa, Giulia Frascaco, Chiara Manganeli, Vincenzo De Santis, and Massimo Pacilli

9	Left and Right Atria: Frequent Imaging in ICU Patients	111
	Luigi Tritapepe, Francesca Pompei, Maurizio Marandola, Giovanni Carriero, Alessandra Di Persio, and Claudio Di Giovanni	
10	Pericardium and Pericardial Disease	117
	Ferdinando Luca Lorini, Stefania Cerutti, and Giovanni Di Dedda	
11	The Aorta: Frequent Imaging in ICU Patients	125
	Luigi Tritapepe, Domenico Vitale, Mario Mezzapesa, Emilia Delia, Ilaria Chiesa, and Roberto Arzilla	
12	The Vena Cava's and the Great Vessels Ultrasound	137
	Silvia Laviola, Angelica Venni, and Silvio Cantini	
13	Ischemia and Myocardial Infarction	145
	Ferdinando Luca Lorini, Marialuigia Dello Russo, and Elena Pagani	
14	The Cardiomyopathies	153
	Ferdinando Luca Lorini, Alessandra Rizza, and Francesco Ferri	
15	Cor Pulmonale and Pulmonary Hypertension	163
	Ferdinando Luca Lorini, Lorenzo Grazioli, and Angelo Vavassori	
16	Mitral Valve	173
	Ilaria Nicoletti, Carla Avallato, and Alessandro Locatelli	
17	Aortic Valve	187
	Alessandro Locatelli, Ilaria Nicoletti, and Carla Avallato	
18	Tricuspid and Pulmonary Valves	201
	Claudio Poli, Armando Sarti, and Vanni Orzalesi	
19	Endocarditis	205
	Roger L. Click	
20	Prosthetic Valve Evaluation	211
	Roger L. Click	
21	Cardiac Tumors and Masses	217
	Roger L. Click	
22	Congenital Septal Abnormalities in The Adult Patient	223
	Ferdinando Luca Lorini, Cristian Ottavio Mirabile, and Moreno Favarato	
23	Essential Pediatric Echocardiography	231
	Ferdinando Luca Lorini, Mariavittoria Lagrotta, and Simona Marcora	

**Part IV Echocardiography in the ICU and OR:
Basic and Advanced Applications**

- 24 Echocardiographic History, Echocardiographic Monitoring, and Goal-Directed, Focus-Oriented, and Comprehensive Examination.** 243
Armando Sarti, Simone Cipani, and Massimo Barattini
- 25 Echocardiography in the ICU and OR: Basic and Advanced Applications** 251
Alessandro Forti, Ferdinando Luca Lorini, and Carlo Sorbara
- 26 General Hemodynamic Assessment.** 257
Carla Avallato, Iliaria Nicoletti, and Alessandro Locatelli
- 27 Contrast Echocardiography in the ICU and OR.** 269
Paolo Voci, Luigi Tritapepe, Marta Iaconi, Daria Collacchi, Demetrio Tallarico, and Luciano Agati
- 28 Echo-Guided Therapy of Myocardial Ischemia.** 275
Michele Oppizzi, Marco Ancona, Vittorio Pazzanese, and Rachele Contri
- 29 Hypovolemia and Fluid Responsiveness.** 283
Armando Sarti, Simone Cipani, and Massimo Barattini
- 30 ARDS, Mechanical Ventilation, and Weaning** 295
Federica Marini, Simone Cipani, Armando Sarti, and Carla Farnesi
- 31 Hypotension.** 305
Luigi Tritapepe, Cecilia Nencini, Giulia Frascacco, and Demetrio Tallarico
- 32 Suspicion of Pulmonary Embolism** 317
Alessandro Locatelli, Iliaria Nicoletti, and Carla Avallato
- 33 Suspicion of Acute Aortic Diseases** 323
Luigi Tritapepe, Francesca Pacini, Giulia Frascacco, Mario Mezzapesa, and Maurizio Caruso
- 34 Chest Pain** 333
Michele Oppizzi, Vittorio Pazzanese, Roberto Spoladore, and Rachele Contri
- 35 Acute Dyspnea.** 351
Simone Cipani and Federica Marini
- 36 Unexplained Hypoxemia** 363
Ferdinando Luca Lorini, Bruno Rossetto, and Francesco Ferri
- 37 Sepsis and Septic Shock** 369
Armando Sarti, Simone Cipani, and Germana Tuccinardi

38	Chest Trauma	375
	Fabio Sangalli, Serena Calcinati, Lucia Galbiati, and Roberto Fumagalli	
39	Acute Atrial Fibrillation and Other Arrhythmias	387
	Piercarlo Ballo	
40	Ultrasound-Guided Nursing	403
	Marco Sozzi, Dorella Donati, and Stella Neri	
41	Multiorgan Donor and Transplanted Patient	413
	Ferdinando Luca Lorini and Lorenzo F. Mantovani	
42	New-Onset Cardiac Murmur in the Unstable Patient	419
	Michele Oppizzi, Marco Ancona, and Vittorio Pazzanese	
43	ICU Echocardiography and Noninvasive Haemodynamic Monitoring	431
	Giulia Frasacco, Mario Mezzapesa, Giovanni Carriero, Fernando Piscioneri, and Luigi Tritapepe	
44	Arterial Ventricular Coupling (AVC)	439
	Simone Cipani, Francesco Maria Traina, and Armando Sarti	
45	Non-cardiac Surgery: Perioperative Echocardiography and Lung Echography	445
	Stefano Romagnoli, Cosimo Chelazzi, and Zaccaria Ricci	
46	Perioperative Echocardiography for Aortic Procedures	453
	Ilaria Blangetti and Alessandro Locatelli	
47	Obstetric Echocardiography	467
	Ferdinando Luca Lorini and Chiara Viviani	
 Part V Ultrasound in the ICU: Other Applications		
48	Echocardiography and Advanced Life Support	477
	Simone Cipani, Rita Cammelli, Sara Felici, and Mauro Cavuta	
49	Ultrasound-Guided Central and Peripheral Vein Cannulation	483
	Fulvio Pinelli and Antonio Franco	
50	Essential Ultrasound for Venous Thrombosis	489
	Federica Marini, Silvia Laviola, Paola Pieraccioni, and Armando Sarti	
51	Lung and Pleural Ultrasonography in Emergency and Intensive Care	495
	Federica Marini, Francesca Covani Frigieri, and Diletta Guarducci	

52	Ultrasound and Morbidly Obese Patient	503
	Pavoni Vittorio, Poli Claudio, and Cipani Simone	
53	Ultrasound Exam Approach in Trauma Patients	513
	Gianfranco Giannasi	
54	Color Duplex Ultrasonography in the Study of the Stenosis of Renal Arteries	533
	Gianfranco Giannasi and Antonio Mannarino	
55	Ultrasound for Percutaneous Tracheostomy	545
	Massimo Barattini, Gabriele Gori, and Alessandra Nella	
56	Transcranial Doppler Ultrasonography in Emergency and Intensive Care	551
	Mauro Pratesi and Daniele Cultrera	
57	Ultrasonography of the Optic Nerve	561
	Daniele Cultrera and Mauro Pratesi	
58	Muscle Ultrasound for the ICU Patient	565
	Duccio Conti and Laura Salucci	
59	Fever of Unknown Origin and Ultrasound	575
	Diletta Guarducci and Armando Sarti	
60	US-Guided Regional Analgesia in the ICU	583
	Stefano Secchi, Juri Valoriani, and Paolo Cecconi	

Part I

Ultrasound and Use of the Echo Machine

Essential Physics of Ultrasound and Use of the Ultrasound Machine

1

Dionisio F. Colella, Paolo Prati, and Armando Sarti

1.1 Ultrasound

Sound is a mechanical wave made up of compressions and rarefactions of molecules in a medium (solid, liquid, or gas) (Fig. 1.1).

Sounds is characterized by some parameters:

- *Frequency* is the number of cycle per unit time (1 s), measured in hertz (Hz). The higher the frequency, the better the resolution, but the lower the penetration (Fig. 1.2).
- *Period* is the duration of a cycle (the inverse of frequency).
- *Wavelength* is the distance that sound travels in one cycle. The wavelength depends on the size of the piezoelectric crystals in the transducer and the medium through which the sound wave travels (Table 1.1).
- *Amplitude* is the amount of change in the oscillating variable. Amplitude decreases as the wave travels (attenuation), leading to echoes from deeper structures being weaker




Fig. 1.1 A sound wave

than those from superficial structures. It is measured in decibels:

$$\text{Decibel}(\text{dB}) = 20 \log_{10} A^2 / A_r^2,$$

where A is the sound amplitude of interest and A_r is a standard reference sound level.

- *Intensity* is the measure of the energy in a sound beam. It is related to potential tissue damage. For example, ultrasound used for lithotripsy has high intensity to fragment renal stones. It is measured in watts per square meter.
- *Power* is the amount of energy transferred. It is expressed in watts.
- The power or the intensity levels are not represented on the ultrasound machine, but there are two other variables that indirectly change those two parameters: mechanical index and

D. F. Colella (✉)
Department of Anesthesia and Intensive Care,
Tor Vergata University, Rome, Italy

P. Prati
Policlinico Tor Vergata, Rome, Italy

A. Sarti
Department of Anesthesia and Intensive Care,
Tor Vergata University, Rome, Italy

Policlinico Tor Vergata, Rome, Italy

Fig. 1.2 Relationship between transducer frequency, penetration, and wavelength. As the transducer frequency increases, resolution increases and penetration decreases




Table 1.1 Relationship between frequency and wavelength

Frequency (MHz)	Wavelength (mm)
1.25	1.2
2.5	0.60
5.0	0.30
7.5	0.20
10.0	0.15

Table 1.2 Ultrasound velocities in different mediums

Material	Velocity (m/s)
Air	330
Water	1497
Fat	1440
Blood	1570
Soft tissue	1540

thermal index. The first one represents the risk of cavitation. The second one is related to the increase of temperature of the tissues (Table 1.1, Fig. 1.2).

- *Propagation velocity* is the velocity determined by the medium that the sound passes through. It is related to the tissue's resistance to compression. Velocity is the product of frequency and wavelength. The propagation velocity through a medium is increased by increasing stiffness of the medium and is reduced by increasing density of the medium (Table 1.2).

Velocity is the product of wavelength and frequency:

$$v = \lambda \times f.$$

1.2 Interaction of Ultrasound with Tissues

1.2.1 Attenuation

When the ultrasound beam passes through uniform tissues, its energy is attenuated by dispersion and absorption.

Absorption is the conversion of ultrasound energy into heat. The attenuation coefficient relates the amount of attenuation to the frequency of the ultrasound beam and the distance that beam travels.

Dispersion occurs because of reflection, refraction, and scattering. The attenuation of the sound wave is increased at higher frequencies, so in order to have better penetration of deeper tissues, a lower frequency is used.

Attenuation involves less energy returning to the transducer, resulting in a poor image.


As the sound traveling through a tissue reaches another tissue with different acoustic properties, the sound energy can be reflected or change its direction, depending on the acoustic impedance of the second interface.

Acoustic impedance is the ability of a tissue to transmit sound and depends on:

- The density of the medium
- The propagation velocity of ultrasound through the medium:

$$Z = \rho \times v,$$

Fig. 1.3 The left ventricular wall is hidden behind a calcified posterior mitral leaflet



where Z is the acoustic impedance, ρ is the density of the material, and v is the speed of ultrasound.

If different mediums have a large difference in acoustic impedance, there is an acoustic impedance mismatch. The greater the acoustic mismatch, the greater the percentage of ultrasound reflected and the lower the percentage transmitted.

1.2.2 Reflection

When a sound wave reaches a smooth surface, it is reflected with an angle that is opposite the incident angle. *The more the angle is near 90° , the lower the amount of energy that is lost.*

There are two types of reflection:

1. *Specular* reflection
2. *Scattering* reflection

If the sound wave reaches a small and irregularly shaped surface (such as red blood cells), the ultrasound energy is scattered in all directions.

Reflection can be measured by the reflection coefficient:

$$R = (Z_2 - Z_1)^2 / (Z_2 + Z_1)^2,$$

where R is the reflection coefficient and Z is the acoustic impedance.

When the second medium encountered is a strong reflector, some phenomena can occur:

- Acoustic shadowing (Fig. 1.3)
- Reverberation (Fig. 1.4)
- A side lobe (Fig. 1.5)

1.2.3 Refraction

When a sound beam reaches the interface between two mediums, some of it is not reflected but passes through the interface, and its direction is altered. This is called refraction. The amount of refraction is proportional to the difference in the velocity of sound in the two tissues and to the angle of incidence:

$$n_1 / n_2 = \sin \theta_1 / \sin \theta_2,$$

where n is the refraction coefficient and θ is the angle of incidence.

It is possible to see some refraction artifact (Figs. 1.6 and 1.7).

Fig. 1.4 Comet tail.
Mirror image: double-
barred aorta





Fig. 1.5 Side lobe artifacts can create a
false aortic flap



1.3 Ultrasound Wave Formation

Ultrasound waves are generated by *piezoelectric crystals*. An electrical current applied to a crystal causes vibration and consequent expansion and contraction. These changes are transmitted into the body as ultrasound waves. Modern transducers are both transmitters and receivers.

There is a strict relationship between time, distance, and velocity of ultrasound propagation.

Knowing the time required for sound to travel from the transducer to an object, the time needed for the returning echo from that object to the transducer, and the propagation velocity in that medium allows one to calculate the distance the ultrasound waves have crossed. This is the basis of ultrasonic imaging.

Electrical energy is not applied to the transducer in a continuous way: ultrasound waves are produced at regular intervals with a pulsed

Fig. 1.6 Grading lobe. The pulmonary catheter seems to be in the aorta





Fig. 1.7 Reflection, refraction, and attenuation

repetition period, leading to a defined *pulse repetition frequency* (PRF; in kilohertz). The wavelength of the ultrasound generated is inversely related to the thickness of the piezoelectric elements.

The piezoelectric elements cannot emit a second pulse until the first has returned to the transducer: the ability to recognize different objects is related to the frequency of emission of the ultrasound wave pulse.

The ultrasound beam emitted from the transducer has a particular shape: it begins with a

narrow beam (near field), and then the ultrasound beam diverges in the far field. The length of the near field (or Fresnel zone) is related to the diameter of the transducer (D) and the wavelength:

$$L_n = D^2 / 4\lambda.$$

Even the angle of divergence, forming the far field (or Fraunhofer zone), is related to the diameter of the transducer (D) and the wavelength:

$$\sin \theta = 1.22\lambda / D.$$

The resolution is improved in the near field because of the narrower diameter of the ultrasound beam. It is easy to understand that a high-diameter transducer with high frequency (short wavelength) can produce the best ultrasound beam.

There is another way to reduce the diameter of the ultrasound beam and thus improve the resolution: focusing the beam. This produces a reduction of the beam size at a particular point, ameliorating the image.

1.3.1 Resolution

This is the ability to recognize two objects. Spatial resolution is the ability to differentiate two separate objects that are close together. Temporal resolution is the ability to place structures at a particular time.

1.3.2 Axial Resolution

This is the ability to recognize two different objects at different depths from the transducer along the axis of the ultrasound beam (Figs. 1.8 and 1.9):




Fig. 1.8 Axial resolution. The spatial pulse length is short enough to be placed within two different structures, so they are resolved




Fig. 1.9 Axial resolution and transducer frequency. Closer objects cannot be resolved by a low transducer frequency. Increasing the transducer frequency (shortening the spatial pulse length and duration) is required to resolve the objects

$$\text{Axial resolution} = \frac{\text{spatial pulse length (SPL)}}{2},$$

where $\text{SPL} = \lambda \times \text{no. of cycles}$.

It is improved by higher-frequency (shorter-wavelength) transducers but at the expense of penetration. Higher frequencies are therefore used to image structures close to the transducer.

1.3.3 Lateral Resolution

This is the ability to distinguish objects that are side by side. It is dependent on the beam width because two objects side by side cannot be distinguished if they are separated by less than the beam width. It is improved by the use of higher-frequency transducer (which increases the beam width) and an optimized focal zone (Figs. 1.10 and 1.11).




Fig. 1.10 Lateral resolution. Wider beams cannot resolve near objects




Fig. 1.11 Lateral resolution. At low depth, lateral resolution is worsened

1.3.4 Temporal Resolution

This is dependent on the frame rate. It is improved by:

- Minimizing depth—the maximum distance from the transducer as this affects the PRF
- Narrowing the sector to the area of interest—narrowing the sector angle
- Minimizing the line density (but at the expense of lateral resolution)

1.4 Doppler Echocardiography

Doppler echocardiography is a method for detecting the direction and velocity of moving blood within the heart.

The *Doppler effect* (or *Doppler shift*) is the change in frequency of a wave for an observer moving relative to the source of the wave (Fig. 1.12).

When the source of the sound wave is moving toward the observer, each successive wave is emitted from a position closer to the observer than the previous wave, and it takes less time than the previous wave to reach the observer. Then the time between the arrival of successive waves is reduced, resulting in a higher frequency. If the source of waves is moving away from the observer, the opposite effect can be seen, with increased time between the arrival of successive waves, giving them a lower frequency.

The amount of that change in frequency is the Doppler shift. Blood flow velocity (V) is related to the Doppler shift by the speed of sound in blood (C) and the intercept angle (θ) between the ultrasound beam and the direction of blood flow:

$$\text{Doppler shift} = 2 \times F(\text{transmitted}) \times [(V \times \cos \theta)] / C$$

A factor of 2 is used because the sound wave has a “round-trip” transit time to and from the transducer. If the ultrasound beam is not parallel to blood flow, an angle of incidence greater than 30° can underestimate the Doppler shift.

There are two kinds of Doppler application: pulsed-wave Doppler and continuous-wave Doppler.

In the *continuous-wave Doppler* technique, the transducer continuously transmits and receives ultrasound waves (Fig. 1.13).

The continuous-wave Doppler technique measures all velocities along the ultrasound beam. It cannot discriminate the time interval from the emission and the reflection, giving no information about the depth of the received signal. The continuous-wave Doppler technique is able to detect very high velocities, and it can be useful to evaluate the high velocity flow through a stenotic aortic valve.

In the *pulsed-wave Doppler* technique, the transducer alternately transmits and receives the ultrasound wave and its returning echo (Fig. 1.14).

Fig. 1.12 Doppler shift




Fig. 1.13 Continuous-wave Doppler imaging





Fig. 1.14 Pulsed-wave Doppler echocardiography





The transducer must wait for the returning echo before sending out another ultrasound wave. The pulsed-wave Doppler technique samples velocities at a specific point (sample volume) of the ultrasound beam.

The number of pulses transmitted from a Doppler transducer each second is called the pulse repetition frequency (PRF). The sampling rate determines the acquisition of information. If the Doppler shift frequency is higher than the PRF, the Doppler signal is displayed on the other side of the baseline. This is the alias artifact.

Aliasing occurs when the measured velocity is greater than half of the PRF (Figs. 1.15 and 1.16). This velocity is called the *Nyquist limit*.

There are some ways to improve the velocity performance of the pulsed-wave Doppler technique:

1. Decrease the distance between the transducer and the sample volume. Reducing the distance the ultrasound beam has to travel will increase the frequency of emission of the pulsed wave (PRF).

Fig. 1.15 Aliasing**Fig. 1.16** Aliasing

2. Choose a low frequency of emission.
3. Set the baseline to display a greater range of velocities (Fig. 1.17).

In *tissue Doppler* imaging (Fig. 1.18), a low-pass filter is used to measure only the velocity of myocardial tissues. Tissue Doppler imaging uses a small pulsed-wave sample volume showing low velocity–high amplitude signals.

Color Doppler imaging combines a 2D image with a Doppler method to visualize the velocity of blood flow within an image plane. The Doppler

shift of thousands of sample volumes displays the directions of the blood cells: blue for away from and red for toward the transducer (Fig. 1.19).

High-flow velocities are displayed in yellow (toward the transducer) and cyan (away from the transducer); green is used to visualize areas of turbulence. As with the pulsed-wave Doppler technique, the color flow Doppler technique suffers from the Nyquist limit, and aliasing can occur (Fig. 1.20).

Color M-mode Doppler imaging combines the spatiotemporal graphic representation of M-mode

Fig. 1.17 Aliasing resolved by setting the right baseline




Fig. 1.18 Tissue Doppler imaging




Fig. 1.19 Color Doppler imaging. *Blue* away from and *red* toward the transducer




Fig. 1.20 Color Doppler aliasing





Fig. 1.21 Color M-mode imaging



and color codification. It shows at the same time a one-dimensional view of anatomic structures and color flow visualization. It is useful to assess transmittal flow (Fig. 1.21).

lar patient to be examined. The positioning of the patient and the probe is discussed in Chap. 2, together with all transthoracic views. First, it is essential to study well the instruction manual of the device at one's disposal to use it optimally.

1.5 Use of the Ultrasound Machine: Optimizing the Picture

The image quality depends on the operator's skill and also on the adjustment of the ultrasound machine according to the features of the particu-

1.5.1 Environment

The brightness of the environment where the examination is done should be reduced. The examination is performed in the ICU at the bedside, so it is preferable to have beds that can be

easily arranged with Trendelenburg positioning, anti-Trendelenburg positioning, head and trunk lifting, and side tilting.

1.5.2 Ultrasonograph Setting

1.5.2.1 Electrocardiogram

Despite often being omitted to save time, it is important to always connect the electrocardiogram (ECG) wire of the ultrasound device to the patient electrodes. The ECG trace, recorded at the base of the display with a “marker” of time coinciding with the moving image, allows establishment of the phases of the cardiac cycle based on electrical activity of the heart, apart from monitoring the ECG. The mechanical systole usually begins immediately after the R wave and ends at about half of the T wave. The end of diastole coincides with the R-wave peak of the ECG (Fig. 1.22).

1.5.2.2 Probes

The probes used for adult echocardiography emit ultrasound with a frequency of about 3 MHz, whereas the probes used in pediatric echocardiography emit higher frequencies, from 5 to 7.5 MHz. These emissions represent the best compromise for use according to different types of patients, given that the higher the frequency, the better the image definition, but the lower the penetration of ultrasound into the tissues. Modern equipment



Fig. 1.23 Ultrasound probes. From left to right: vascular and soft tissue linear probe, cardiac phased-array probe, abdominal convex probe

can produce sharper images through “tissue harmonic imaging”: in a nutshell, the harmonics produced by the interaction of the ultrasound beam with the tissues are enhanced, and the fundamental harmonic frequencies are suppressed, resulting in better far-field quality. The image quality is better, but the very echo-reflective structures, such as pericardium and valves, can thus appear thicker than they really are. The probes have a touch and often light marker, defining the scanning plane and laterality. Figure 1.23 shows the various probes used for the study of the heart with trans-thoracic approaches, the vessels, and the thoracic and abdominal organs.

1.5.2.3 Sector Depth

The depth can be adjusted by the operator. The machine starts with a default standard depth so that the whole heart is displayed, but the depth of the field can be varied in order to position the structures of interest in the middle of the image. If the outer edges of the heart in the default image exceed the limits of the display, the heart as a whole or some part of it is certainly enlarged.

1.5.2.4 Width of the Scanning Beam

The maximum amplitude ensures that the most lateral structures are seen, but a reduction of the amplitude may sometimes be preferable to produce greater definition of the central structures; this is because a shorter time is needed to scan a narrower angle.




Fig. 1.22 ECG systolic and diastolic phases

1.5.2.5 Gain

The construction of the image as grayscale or monochrome images depends on the intensity of the return signal, which depends on the distance traveled and on the reflective properties of the tissues encountered (see earlier). Therefore, the gain can be adjusted for different depths in order to compensate for the reduction of the return signal. This adjustment (time gain compensation), which is automatic in modern equipment, usually occurs through a system of levers that correspond to vertical depths of the field. Observing the display as the default, the operator improves the image manually by moving the levers that correspond to different levels of depth. In some devices there is also a system of horizontal adjustment for adjusting the image in the lateral fields (lateral gain compensation). However, these adjustments must be done with care since excess gain produces brighter images, leading to poor definition between close structures and even artifacts. In contrast, too dark images are produced with not enough gain, hiding some low-echo-reflective structures. Even though it depends somewhat on operator preference, a well-adjusted image (Fig. 1.24) is one that has:

- Fairly uniform intensity of solid structures
- A slight speckling in the dark cavities full of blood

1.5.2.6 Focus

The focus of the image is usually by default in the central part of the display, but one can move the focus to higher or lower levels for further research on particular structures.

1.5.2.7 Regulating Continuous-Wave Doppler and Pulsed-Wave Doppler Imaging

As already mentioned, continuous-wave Doppler imaging is used for the measurement of high flow rate in line with the cursor all along the stream to be examined. Pulsed-wave Doppler imaging is not suitable for high-speed flows but reproduces the flow in a specific area to be examined. The operator can adjust the gain of the Doppler signal. The optimal image is one that shows well the shape of the wave flow (changes in speed over time). By convention, the blood flow movement toward the transducer is represented above the baseline. In contrast, the movement away from the transducer is represented below the baseline. The scale of reproduction of the Doppler signal (y -axis) can be adjusted to avoid cutting high-speed-flow waves. The speed, usually 50 mm/s (x -axis), can be adjusted to better fit the times for special measures. The alignment of the ultrasound beam with the flow remains paramount for both continuous-wave and pulsed-wave interrogation. An angle of more than 30° between the blood flow

Fig. 1.24 Echocardiographic image (apical five-chamber view) with good contrast adjustment

